

Clear identification, cleanliness, sterilizability, and permanence are offered with AdvantaPure's AdvantaLABEL tagging system. AdvantaLABEL eliminates the issues involved with old-style metal, plastic, or even paper tags. Since the label is sealed to the hose in clean silicone, it can be either steam sterilized (SIP) or chemically cleaned (CIP) many times without degradation of the label.



AdvantaLABEL is custom configured to your requirements – part number, manufactured date, planned replacement date, or any specific information you deem critical. Contact your AdvantaPure representative today to see how AdvantaLABEL can solve your hose tagging needs.



## KEY FEATURES

- O.D. mounted to reinforced silicone hoses – does not contact material flow
- Permanently vulcanized silicone
- Certified free of animal-derived ingredients
- Easy to read – remains bright after repeated sterilization
- Autoclavable
- Zero entrapment
- Handles clean-in-place (CIP) or steam-in-place (SIP) processes
- Customer supplied text is application specific
- Add color striping for quick line identification – 10 standard colors available\*
- Field installation methods available



IDENTIFICATION

www.advantapure.com

- RED
- ORANGE
- YELLOW
- GREEN
- BLUE
- PURPLE
- BROWN
- BLACK
- GRAY
- WHITE

\* Actual colors may vary



Customize AdvantaLABEL to include the specific data you need.

# MORE ADVANTAPURE PRODUCTS

- Unreinforced platinum silicone tubing
- High pressure platinum silicone tubing
- Ultra low temperature platinum silicone tubing
- Braid-reinforced platinum silicone hose
- Stainless steel wire and polyester mesh reinforced wrapped platinum silicone hose
- AdvantaFlex® BioPharmaceutical Grade Tubing
- PTFE hose with stainless steel overbraiding (with or without overjacketing)
- Rubber-covered EPDM hose
- Stainless steel sanitary crimp fittings and complete hose assemblies including molded silicone ends
- Custom hose identification solutions – laser-etched stainless steel collars, tubing/hose printing, and solid color exterior layers
- Identification solutions including GammaTag®, Hose Track® and P•E•T Process Equipment Tracking® from NewAge® Industries' Verigenics® division
- Special hose coverings, cleaning, and packaging



## Purity in Fluid Flow Systems®



Made in U.S.A.

Address: NewAge Industries AdvantaPure®  
145 James Way  
Southampton, PA 18966 U.S.A.



Toll Free Phone: **888-755-4370**  
Phone: **215-526-2151**  
Toll Free Fax: **888-258-4293**  
Fax: **215-526-2167**  
Web Site: **www.advantapure.com**  
E-mail: **sales@advantapure.com**

DISTRIBUTED BY:

©2018 NewAge® Industries, Inc. Specifications are subject to change without notice. AdvantaPure products are produced by NewAge® Industries, Inc. AdvantaLabel®, AdvantaFlex®, GammaTag®, Hose Track®, P•E•T Process Equipment Tracking®, Purity in Fluid Flow Systems®, NewAge Industries AdvantaPure®, NewAge®, and the "N" logo are trademarks of NewAge® Industries, Inc.



WARRANTY: NewAge Industries AdvantaPure products are warranted to be free from defects in material and workmanship. Any product found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of NewAge Industries AdvantaPure products, item for item, which after examination by NewAge Industries AdvantaPure are deemed defective. NewAge Industries AdvantaPure specifically disclaims any other liability. The data supplied is provided as a helpful guideline and is believed to be reliable; however, nothing stated shall constitute a guarantee, recommendation, or warranty for any application.  
NOTICE: NewAge Industries, Inc. has been a privately owned corporation since 1954. Neither the corporation nor its principals are affiliated, in any way, with any other group or entity either representing themselves as or being referred to by the term New-Age.