

Clear identification, cleanliness, sterilizability, and permanence are offered with AdvantaPure's AdvantaLABEL tagging system. AdvantaLABEL eliminates the issues involved with old-style metal, plastic, or even paper tags. Since the label is sealed to the hose in clean silicone, it can be either steam sterilized (SIP) or chemically cleaned (CIP) many times without degradation of the label.



AdvantaLABEL is custom configured to meet your specific needs – part number, manufactured date, planned replacement date, or any specific information you deem critical. Contact your AdvantaPure representative today to see how AdvantaLABEL can solve your hose tagging needs.



KEY FEATURES




- O.D. mounted to reinforced silicone hoses – does not contact material flow
- Permanently vulcanized silicone
- Easy to read – remains bright after repeated sterilization
- Autoclavable
- Zero entrapment
- Handles clean-in-place (CIP) or steam-in-place (SIP) processes
- Customer supplied text is application specific
- Add color striping for quick line identification – 10 standard colors available*
- Field installation methods available

IDENTIFICATION

www.advantapure.com

- RED
- ORANGE
- YELLOW
- GREEN
- BLUE
- PURPLE
- BROWN
- BLACK
- GRAY
- WHITE



* Actual colors may vary

Customize AdvantaLABEL to include the specific data you need most.

What Products Does AdvantaPure Offer?

AdvantaPure's prime product offerings consist of low-volatile platinum silicone tubing and hose meeting USP Class VI, FDA CFR 177.2600, and ISO 10993, and overbraided PTFE hose.

Our line of products includes:

- Unreinforced 50 Shore A platinum silicone tubing
- Polyester braid-reinforced platinum silicone hose
- Stainless-steel wire and polyester mesh-reinforced wrapped platinum silicone hose
- Polyester mesh-reinforced 4-ply platinum silicone hose
- PTFE hose with stainless-steel or polypropylene overbraiding (with or without armor casing or silicone jacketing)
- Rubber-covered FEP or EPDM hose
- Stainless-steel sanitary crimp fittings and complete hose assemblies
- Various custom hose identification solutions including silicone-encased labels with color striping, laser-etched stainless steel collars, tubing/hose printing, and colored and color-striped silicone tubing, colored braid material for hose, and RFID tags.
- Special cleaning and packaging



Purity in Fluid Flow SystemsSM



Address: NewAge® Industries AdvantaPure®
145 James Way
Southampton, PA 18966 U.S.A.

Toll Free Phone: **888-755-4370**
Phone: **215-526-2151**
Toll Free Fax: **888-258-4293**
Fax: **215-526-2167**
Web Site: **www.advantapure.com**
E-mail: **sales@advantapure.com**

AdvantaPURE
NEW AGE INDUSTRIES



DISTRIBUTED BY:

©2008 NewAge® Industries, Inc.
AdvantaPure® products are produced by NewAge® Industries, Inc.
AdvantaPure® and the AdvantaPure® logo are registered trademarks of NewAge® Industries, Inc.
NewAge® and the "N" logo are registered trademarks of NewAge® Industries, Inc.
Purity in Fluid Flow Systems™ is a service mark of NewAge® Industries, Inc.



ESOP Employee Owned
for Your Benefit

Celebrating over
50
Years
of Service
1954-2008

WARRANTY: NewAge Industries AdvantaPure products are warranted to be free from defects in material and workmanship. Any product found to carry such defects will be replaced without charge. This warranty is limited only to the replacement of NewAge Industries AdvantaPure products, item for item, which after examination by NewAge Industries AdvantaPure are deemed defective. NewAge Industries AdvantaPure specifically disclaims any other liability. The data supplied is provided as a helpful guideline and is believed to be reliable; however, nothing stated shall constitute a guarantee, recommendation, or warranty for any application.

NOTICE: NewAge Industries, Inc. has been a privately owned corporation since 1954. Neither the corporation nor its principals are affiliated, in any way, with any other group or entity either representing themselves as or being referred to by the term New-Age.