

How to Order an AdvantaPure Hose Assembly

Example:



This assembly consists of a 3/4" I.D. double-braided silicone hose with a 3/4" mini sanitary Tri-Clamp® fitting (316L stainless steel) on one end and a 1-1/2" sanitary Tri-Clamp step up fitting (316L stainless steel) on the other. The assembly is 72" in length and uses an AdvantaLabel® for hose identification.

Collars are used to attach fittings to hoses. When pricing an assembly, remember to add the cost for all collars.

Hose I.D. Size

Code	Size
02	1/8"
04	1/4"
06	3/8"
08	1/2"
10	5/8"
12	3/4"
16	1"
24	1-1/2"
32	2"
40	2-1/2"
48	3"
64	4"

Hose Style

Code
APAF-BP
APST
APSPG
APSH-P
APSH-DB
APSW-P
APSW-PC
APFOS-W
APFOS-WC
APFOSJ-W
APFOSJ-WC
APFOP-WC
APFRC-W-FLX
APWF-W

Fitting Type

Code	Description
01	Male Pipe
015	Male Pipe Step Down
02	Female Pipe
025	Female Pipe Step Down
07	Female JIC Swivel
10	Tri-Clamp (Sanitary)
101	Tri-Clamp (Sanitary) Step Up
102	Tri-Clamp (Sanitary) 2X Step Up
11	Mini Tri-Clamp (Sanitary)
111	Mini Tri-Clamp (Sanitary) Step Up
115	Mini Tri-Clamp (Sanitary) Step Down
13	Tri-Clamp, Flare Thru
16	Flange Retainer w/150 lb. Flange
17	Flange Retainer w/300 lb. Flange
18	Male Cam & Groove
22	Female Cam & Groove Swivel (locking arms only)
25	Female Cam & Groove Swivel, Flare Thru
26	Male I-Line
27	Female I-Line
28	Flange Retainer, Flare Thru with 150 lb. Flange
30	Male Bevel Seat
31	Female Bevel Seat Assembly with Nut
37	Tube Weld Adapter
38	Pipe Weld Adapter
40	Compression Tube Adapter
42	90 Degree Mini Tri-Clamp (Sanitary)
421	90 Degree Mini Tri-Clamp (Sanitary) Step Up
43	45 Degree Mini Tri-Clamp (Sanitary)
46	90 Degree Tri-Clamp (Sanitary)
461	90 Degree Tri-Clamp (Sanitary) Step Up
47	45 Degree Tri-Clamp (Sanitary)
53	Male Swivel Union Assembly
54	Female Swivel Union Assembly

Fitting Material

Code	Description
B	Molded with 304 Stainless Steel Backing Cup
C	Molded with Polysulfone Backing Cup
E	Powder Coated Carbon Steel
F	Fluoropolymer Encapsulated SS
H	Hastelloy
M	Monel
N	AL6XN
R	304 Stainless Steel
S	316L Stainless Steel
W	Passivated Stainless Steel
X	Electropolished Stainless Steel
Y	Overmolded Polypropylene
Z	PVDF

Accessories & Services

Code	Description
75	Armor Casing, 300 Series Stainless Steel
76	Armor Cuffs, 300 Series Stainless Steel (specify length)
77	Firesleeve
79	Heat Shrink Polyolefin Cover
80	Stainless Steel Overbraid
81	Silicone Cover
82	AdvantaLabel - Silicone
83	Laser Etched Collars
84	Laser Etched Tags
GS	Gamma Sterilized

3/17/17



PURITY IN FLUID FLOW SYSTEMS®

888-755-4370 • 215-526-2151

www.advantapure.com