

# How to Order an AdvantaPure Hose Assembly

## Example:



This assembly consists of a 3/4" I.D. double-braided silicone hose with a 3/4" mini sanitary Tri-Clamp® fitting (316L stainless steel) on one end and a 1-1/2" sanitary Tri-Clamp step up fitting (316L stainless steel) on the other. The assembly is 72" in length and uses an AdvantaLabel® for hose identification.

Collars are used to attach fittings to hoses. When pricing an assembly, remember to add the cost for all collars.

## Hose I.D. Size

| Code | Size   |
|------|--------|
| 02   | 1/8"   |
| 04   | 1/4"   |
| 06   | 3/8"   |
| 08   | 1/2"   |
| 10   | 5/8"   |
| 12   | 3/4"   |
| 16   | 1"     |
| 24   | 1-1/2" |
| 32   | 2"     |
| 40   | 2-1/2" |
| 48   | 3"     |
| 64   | 4"     |

## Hose Style

| Code      |
|-----------|
| APAF-BP   |
| APST      |
| APHP      |
| APSPG     |
| APSH-P    |
| APSH-DB   |
| APSW-P    |
| APSW-PC   |
| APSW-SP   |
| APFOS-W   |
| APFOS-WC  |
| APFOSJ-W  |
| APFOSJ-WC |
| APEWF-W   |

## Fitting Type

| Code | Description                                    |
|------|--|
| 01   | Male Pipe                                      |
| 015  | Male Pipe Step Down                            |
| 02   | Female Pipe                                    |
| 025  | Female Pipe Step Down                          |
| 07   | Female JIC Swivel                              |
| 10   | Tri-Clamp (Sanitary)                           |
| 101  | Tri-Clamp (Sanitary) Step Up                   |
| 102  | Tri-Clamp (Sanitary) 2X Step Up                |
| 11   | Mini Tri-Clamp (Sanitary)                      |
| 111  | Mini Tri-Clamp (Sanitary) Step Up              |
| 115  | Mini Tri-Clamp (Sanitary) Step Down            |
| 16   | Flange Retainer w/150 lb. Flange               |
| 17   | Flange Retainer w/300 lb. Flange               |
| 18   | Male Cam & Groove                              |
| 22   | Female Cam & Groove Swivel (locking arms only) |
| 26   | Male I-Line                                    |
| 27   | Female I-Line                                  |
| 30   | Male Bevel Seat                                |
| 31   | Female Bevel Seat Assembly with Nut            |
| 37   | Tube Weld Adapter                              |
| 38   | Pipe Weld Adapter                              |
| 40   | Compression Tube Adapter                       |
| 42   | 90 Degree Mini Tri-Clamp (Sanitary)            |
| 421  | 90 Degree Mini Tri-Clamp (Sanitary) Step Up    |
| 43   | 45 Degree Mini Tri-Clamp (Sanitary)            |
| 46   | 90 Degree Tri-Clamp (Sanitary)                 |
| 461  | 90 Degree Tri-Clamp (Sanitary) Step Up         |
| 47   | 45 Degree Tri-Clamp (Sanitary)                 |
| 53   | Male Swivel Union Assembly                     |
| 54   | Female Swivel Union Assembly                   |

## Fitting Material

| Code | Description                                 |
|------|---|
| B    | Molded with 304 Stainless Steel Backing Cup |
| C    | Molded with Polysulfone Backing Cup         |
| E    | Powder Coated Carbon Steel                  |
| F    | Fluoropolymer Encapsulated SS               |
| H    | Hastelloy                                   |
| M    | Monel                                       |
| N    | AL6XN                                       |
| R    | 304 Stainless Steel                         |
| S    | 316L Stainless Steel                        |
| W    | Passivated Stainless Steel                  |
| X    | Electropolished Stainless Steel             |
| Y    | Overmolded Polypropylene                    |
| Z    | PVDF  |

## Accessories & Services

| Code | Description  |
|------|--|
| 75   | Armor Casing, 300 Series Stainless Steel                 |
| 76   | Armor Cuffs, 300 Series Stainless Steel (specify length) |
| 77   | Firesleeve   |
| 79   | Heat Shrink Polyolefin Cover                             |
| 80   | Stainless Steel Overbraid                                |
| 82   | AdvantaLabel - Silicone                                  |
| 83   | Laser Etched Collars                                     |
| 84   | Laser Etched Tags  |
| GS   | Gamma Sterilized   |

4/9/18



PURITY IN FLUID FLOW SYSTEMS®

888-755-4370 • 215-526-2151

www.advantapure.com