

# How to Order an AdvantaPure Hose Assembly

## Example:



This assembly consists of a 3/4" I.D. double-braided silicone hose with a 3/4" mini sanitary Tri-Clamp® fitting (316L stainless steel) on one end and a 1-1/2" sanitary Tri-Clamp step up fitting (316L stainless steel) on the other. The assembly is 72" in length and uses an AdvantaLabel® for hose identification.

Collars are used to attach fittings to hoses. When pricing an assembly, remember to add the cost for all collars.

## Hose I.D. Size

Code	Size
02	1/8"
04	1/4"
06	3/8"
08	1/2"
10	5/8"
12	3/4"
16	1"
24	1-1/2"
32	2"
40	2-1/2"
48	3"
64	4"

## Hose Style

Code
APSH-P
APSH-DB
APSW-SP
APSW-SPC
APFOS-W
APFOS-WC
APFOSJ-W
APFOSJ-WC
APEWF-W
GSPT

## Fitting Type

Code	Description
01	Male Pipe
015	Male Pipe Step Down
02	Female Pipe
025	Female Pipe Step Down
07	Female JIC Swivel
10	Tri-Clamp (Sanitary)
101	Tri-Clamp (Sanitary) Step Up
102	Tri-Clamp (Sanitary) 2X Step Up
11	Mini Tri-Clamp (Sanitary)
111	Mini Tri-Clamp (Sanitary) Step Up
115	Mini Tri-Clamp (Sanitary) Step Down
16	Flange Retainer w/150 lb. Flange
17	Flange Retainer w/300 lb. Flange
18	Male Cam & Groove
22	Female Cam & Groove Swivel (locking arms only)
26	Male I-Line
27	Female I-Line
30	Male Bevel Seat
31	Female Bevel Seat Assembly with Nut
34	Hose Mender
37	Tube Weld Adapter
38	Pipe Weld Adapter
40	Compression Tube Adapter
42	90 Degree Mini Tri-Clamp (Sanitary)
421	90 Degree Mini Tri-Clamp (Sanitary) Step Up
43	45 Degree Mini Tri-Clamp (Sanitary)
46	90 Degree Tri-Clamp (Sanitary)
461	90 Degree Tri-Clamp (Sanitary) Step Up
47	45 Degree Tri-Clamp (Sanitary)
53	Male Swivel Union Assembly
54	Female Swivel Union Assembly

## Fitting Material

Code	Description
B	Molded with 304 Stainless Steel Backup Cup
C	Molded with Polysulfone Backup Cup
H	Hastelloy
M	Monel
N	AL6XN
R	304 Stainless Steel
S	316L Stainless Steel
X	Electropolished Stainless Steel
Y	Overmolded Polypropylene

## Accessories & Services

Code	Description
75	Armor Casing, 300 Series Stainless Steel
77	Firesleeve
79	Heat Shrink Polyolefin Cover
82	AdvantaLabel - Silicone
83	Laser Etched Collars
84	Laser Etched Tags
GI	Gamma Irradiated

05/05/22



PURITY IN FLUID FLOW SYSTEMS®

888-755-4370 • 215-526-2151

www.advantapure.com